

OPERATION OF THERMOSCIENTIFIC EUTECH COND 6+.

1. Equipment is switched on and let to warm for 10 minutes.
2. Conductivity cell is immersed in 0.01N KCl.
3. The equipment is calibrated by pressing the 'Cal/Back' button.
4. The equipment will display the letters 'CA'.
5. Once the calibration is successful, it will automatically start blinking 1413 μ S.
6. Calibration is completed by pressing the 'Enter/Hold' button.
7. Rinse the electrode with distilled water before starting the analysis.
8. Rinse the electrode with distilled water between samples to avoid contamination and thoroughly after measurements are completed.
9. Switch off the equipment.
10. Dry storage of the electrode is recommended.
11. Any spillages on the meter and its surroundings are cleaned immediately after use.

NB:

Do NOT remove probe guard during measurements.